

Product datasheet for TP301496M

C14orf119 (NM_017924) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 14 open reading frame 119 (C14orf119), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201496 protein sequence Red =Cloning site Green =Tags(s) MPLESSSSMPLSFPSLLPSVPHNTNPSPPLMSYITSQEMKCILHWFANWSGPQRRERFLEDLVAKAVPEKL QPLLDLSLEQLSVSGADRPPSIFECQLHLWDQWFRGWAEQERNEFVRQLEFSEPDFVAKFYQAVAATAGK D TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_060394</u>
Locus ID:	55017
UniProt ID:	<u>Q9NWX9</u>
RefSeq Size:	2960


[View online »](#)

Cytogenetics: 14q11.2

RefSeq ORF: 420

Product images:



Coomassie blue staining of purified C14orf119 protein (Cat# [TP301496]). The protein was produced from HEK293T cells transfected with C14orf119 cDNA clone (Cat# [RC201496]) using MegaTran 2.0 (Cat# [TT210002]).