

Product datasheet for **TP301507M**

ACAD8 (NM_014384) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human acyl-Coenzyme A dehydrogenase family, member 8 (ACAD8), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC201507 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MLWSGCRRFGARLGCLPGGLRVLVQTGHRSLTSCIDPSMGLNEEQKEFQKVAFDFAAREMAPNMAEWDQK
ELFPVDVMRKAQQLGFGVYIQTDVGGSLRLDTSVIFEALATGCTSTTAYISIHNMCAWMIDSFGNEE
QRHKFCPLCTMEKFASYCLTEPGSGDAASLLTSAKKQGDHYILNGSKAFISGAGESDIYVVMCRTGGL
GPKGISCIVVEKGPGLSFGKKEKKVGNWSQPTRAVIFEDCAVPVANRIGSEGQGFLIAVRGLNGGRINI
ASCSLGAHASVILTRDHLNVRKQFGEPLASNQYLQFTLADMATRLVAARLMVRNAAVALQEERKDAVAL
CSMAKLFATDECFaicnQALQMhGGYgYLkDYAVQYVRDSRVHQILEGSNEVMRILISRSLQEQE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	44.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055199



[View online »](#)

Locus ID: 27034

UniProt ID: [Q9UKU7](#)

RefSeq Size: 2216

Cytogenetics: 11q25

RefSeq ORF: 1245


Synonyms: ACAD-8; ARC42; IBDH

Summary: This gene encodes a member of the acyl-CoA dehydrogenase family of enzymes that catalyze the dehydrogenation of acyl-CoA derivatives in the metabolism of fatty acids or branch chained amino acids. The encoded protein is a mitochondrial enzyme that functions in catabolism of the branched-chain amino acid valine. Defects in this gene are the cause of isobutyryl-CoA dehydrogenase deficiency.[provided by RefSeq, Nov 2009]

Protein Families: Transcription Factors

Protein Pathways: Metabolic pathways, Valine, leucine and isoleucine degradation

Product images:



Coomassie blue staining of purified ACAD8 protein (Cat# [TP301507]). The protein was produced from HEK293T cells transfected with ACAD8 cDNA clone (Cat# [RC201507]) using MegaTran 2.0 (Cat# [TT210002]).