

Product datasheet for **TP301557L**

SSNA1 (NM_003731) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Sjogren syndrome nuclear autoantigen 1 (SSNA1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201557 protein sequence Red =Cloning site Green =Tags(s)
	MTQQGAALQNYNNELVKCIEELCQKREELCRQIQEEDEKQRLQNEVRQLTEKLARVNENLARKIASRNE FDRTIAETEAAAYLKILESSQTLLSVLKREAGNLTKATAPDQKSSGGRDS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003722
Locus ID:	8636
UniProt ID:	O43805 , A0A024R8G6
RefSeq Size:	954
Cytogenetics:	9q34.3




[View online »](#)

RefSeq ORF: 357

Synonyms: N14; NA-14; NA14

Product images:



Coomassie blue staining of purified SSNA1 protein (Cat# [TP301557]). The protein was produced from HEK293T cells transfected with SSNA1 cDNA clone (Cat# [RC201557]) using MegaTran 2.0 (Cat# [TT210002]).