

Product datasheet for TP301561

GALE (NM_000403) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human UDP-galactose-4-epimerase (GALE), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201561 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAEKVLVTGGAGYIGSHTVLELLEAGYLPVIDNFHNAFRGGGSLPESLRRVQELTGRSVEFEEMDILDQ GALQRLFKKYSFMAVIHFAGLKAVGESVQKPLDYRVNLTGTIQLLEIMKAHGKLVFSSATVYGNPQ YLPLDEAHPTGGCTNPYGKSKFFIEEMIRDLQCADKTWNAVLLRYFNPTGAHASGCIGEDPQGIPNNLMP YVSQVAIGRREALNVFGNDYDTEGTGVRDYIHVVDLAKGHIAALRKLKEQCGCRIYNLGTGTGYSVLQM VQAMEKASGKKIPYKVARREGDVAACYANPSLAQEELGWTAALGLDRMCEDLWRWQKQNPSTSGFGTQA</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	38.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000394
Locus ID:	2582



[View online »](#)

UniProt ID: [Q14376](#), [A0A384NL38](#)

RefSeq Size: 1647

Cytogenetics: 1p36.11

RefSeq ORF: 1044


Synonyms: SDR1E1

Summary: This gene encodes UDP-galactose-4-epimerase which catalyzes two distinct but analogous reactions: the epimerization of UDP-glucose to UDP-galactose, and the epimerization of UDP-N-acetylglucosamine to UDP-N-acetylgalactosamine. The bifunctional nature of the enzyme has the important metabolic consequence that mutant cells (or individuals) are dependent not only on exogenous galactose, but also on exogenous N-acetylgalactosamine as a necessary precursor for the synthesis of glycoproteins and glycolipids. Mutations in this gene result in epimerase-deficiency galactosemia, also referred to as galactosemia type 3, a disease characterized by liver damage, early-onset cataracts, deafness and cognitive disability, with symptoms ranging from mild ('peripheral' form) to severe ('generalized' form). Multiple alternatively spliced transcripts encoding the same protein have been identified. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Amino sugar and nucleotide sugar metabolism, Galactose metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified GALE protein (Cat# TP301561). The protein was produced from HEK293T cells transfected with GALE cDNA clone (Cat# [RC201561]) using MegaTran 2.0 (Cat# [TT210002]).