

Product datasheet for **TP301563**

EIF2S2 (NM_003908) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human eukaryotic translation initiation factor 2, subunit 2 beta, 38kDa (EIF2S2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201563 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSGDEMIFDPTMSKKKKKKKPFMLDEEGDTQTEETQPSETKEVEPEPTEDKDLEADEEDTRKKDASDDL DDLNFFNQKKKKKKTKKIFDIDEAEEGVKDLKIESDVQEPTPEDDLDIMLGNNKKKKKNVKKPDEDEIL EKDEALEDEDNKKDDGISFSNQTGPAWAGSERDYTYDELLNRVFNIMREKNPDMVAGEKRKFVMPKPQV V RVGTTKTSFVNFTDICKLLHRQPKHLLAFLAELGTSGSIDGNNQLVIKGRFQQKQIENVLRRYIKEYVT CHTCRSPDTILQKDTRLYFLQCETCHSRCSVASIKTGFQAVTGKRAQLRAKAN</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	38.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_003899</u>



[View online »](#)

Locus ID:	8894
UniProt ID:	<u>P20042</u>
RefSeq Size:	2592
Cytogenetics:	20q11.22
RefSeq ORF:	999
Synonyms:	eIF-2-beta; EIF2; EIF2B; EIF2beta; PPP1R67
Summary:	Eukaryotic translation initiation factor 2 (EIF-2) functions in the early steps of protein synthesis by forming a ternary complex with GTP and initiator tRNA and binding to a 40S ribosomal subunit. EIF-2 is composed of three subunits, alpha, beta, and gamma, with the protein encoded by this gene representing the beta subunit. The beta subunit catalyzes the exchange of GDP for GTP, which recycles the EIF-2 complex for another round of initiation. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015]

Product images:



Coomassie blue staining of purified EIF2S2 protein (Cat# TP301563). The protein was produced from HEK293T cells transfected with EIF2S2 cDNA clone (Cat# [RC201563]) using MegaTran 2.0 (Cat# [TT210002]).