

Product datasheet for **TP301590L**

NFIB (NM_005596) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nuclear factor I/B (NFIB), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201590 protein sequence Red=Cloning site Green=Tags(s)

MMYSPICLTQDEFHPFIEALLPHVRAIAYTWFNLQARKRKYFKKHEKRMSKDEERAVKDELLSEKPEIKQ
 KWASRLAKLRKDIRQEYREDFVLTVTGKKHPCCVLSNPDQKGKIRRIDCLRQADKVVRLDLVMVILFKG
 IPLESTDGERLMKSPHCTNPALCVQPHHITVSVKELDLFLAYVQEQDSGQSGSPSHNDPAKNPPGYLED
 SFVKSGVFNVSELVRVSRTPITQGTGVNFPIGEIPSQPYHDMMSGVNLQRSLSPPSSKRPKTISIDEN
 MEPSPTGDFYPSPPSPAAGSRTWHERDQDMSSPTTMKKPEKPLFSSASPQDSSPRLSTFPQHHPGIPGV
 AHSVISTRTPPPSPLPFPTQAILPPAPSSYFSHTIRYPPHLNPQDTLKNYVPSYDPSSPQTSQSWYLG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	47.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005587
Locus ID:	4781



[View online »](#)

UniProt ID:	<u>O00712, Q5VW28</u>
RefSeq Size:	8285
Cytogenetics:	9p23-p22.3
RefSeq ORF:	1260
Synonyms:	CTF; HMGIC/NFIB; MACID; NF-I/B; NF1-B; NFI-B; NFI-RED; NFIB2; NFIB3
Summary:	Recognizes and binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified NFIB protein (Cat# [TP301590]). The protein was produced from HEK293T cells transfected with NFIB cDNA clone (Cat# [RC201590]) using MegaTran 2.0 (Cat# [TT210002]).