

## Product datasheet for TP301620

### DNAJB6 (NM\_005494) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Dnaj (Hsp40) homolog, subfamily B, member 6 (DNAJB6), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201620 protein sequence Red=Cloning site Green=Tags(s)
	MVDYYEVLGVQRHASPEDIKKAYRKLALKWHPDKNPENKEEAERKFKQVAEAYEVLSDAKKRDIYDKYGK EGLNGGGGGGSHFDSPFEGFTFRNPDDVFREFFGGRPFSFDFFEDPFEDFFGNRRGPRGSRSRGTGSF FSAFSGFPSFGSGFSSFDTGFTSFGSLGHGGLTSFSSTSGGSGMGNFKSISTSTKMNVRKITTKRIVE NGQERVEVEEDGQLKSLTINGKEQLLRLDNK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	26.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_005485</a>
Locus ID:	10049



[View online »](#)

UniProt ID: [O75190](#)

RefSeq Size: 1568


Cytogenetics: 7q36.3

RefSeq ORF: 723

Synonyms: DJ4; Dnaj; HHDJ1; HSJ-2; HSJ2; LGMD1D; LGMD1E; LGMDD1; MRJ; MSJ-1

**Summary:** This gene encodes a member of the DNAJ protein family. DNAJ family members are characterized by a highly conserved amino acid stretch called the 'J-domain' and function as one of the two major classes of molecular chaperones involved in a wide range of cellular events, such as protein folding and oligomeric protein complex assembly. This family member may also play a role in polyglutamine aggregation in specific neurons. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been fully described. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified DNAJB6 protein (Cat# TP301620). The protein was produced from HEK293T cells transfected with DNAJB6 cDNA clone (Cat# [RC201620]) using MegaTran 2.0 (Cat# [TT210002]).