

## Product datasheet for TP301680M


### C20orf195 (FNDC11) (NM\_024059) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 20 open reading frame 195 (C20orf195), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201680 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MSTHVAGLGLDKMKLGNPQSFLDQEEADDQQLLEPEAWKTYTERRNALREFLTSDLSPHLLKRHHARM QL LRKCSYYIEVLPHKHLALGDQNPLVLPSALFQLIDPWKFQRMKKVGTAQTKIQLLLGLLEQLDHGRAEL DALLRSPDPRPFLADWALVERRLADVSAVMDSFLTMMVPGRLHVKHRLVSDVSATKIPHIWMLSTKMP V VFDRKASAAHQDWARLRWFVTIQPATSEQYELRFRLLDPRTQQECAQCGVIPVAACTFDVRNLLPNRSYK FTIKRAETSTLVYEPWRDSLTLHTKPEPLEGPALSHSV  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span>
Tag:	C-Myc/DDK
Predicted MW:	36.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_076964</a></u>


[View online »](#)

Locus ID: 79025  
 UniProt ID: [Q9BVV2](#)  
 RefSeq Size: 1124  
 Cytogenetics: 20q13.33  
 RefSeq ORF: 954  
 Synonyms: C20orf195

**Product images:**


Coomassie blue staining of purified FNDC11 protein (Cat# [TP301680]). The protein was produced from HEK293T cells transfected with FNDC11 cDNA clone (Cat# [RC201680]) using MegaTran 2.0 (Cat# [TT210002]).