

## Product datasheet for TP301710M

### PDIA6 (NM\_005742) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein disulfide isomerase family A, member 6 (PDIA6), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201710 protein sequence Red=Cloning site Green=Tags(s)

MALLVLGLVSCTFFLAVNGLYSSSDVIELTPSNFNREVIQSDSLWLVEFYAPWCGHCQRLTPEWKKAAT  
ALKDWKVGAVDADKHHSLGGQYGVQGFPTIKIFGNSKNRPEDYQGGRTGEAIVDAALSALRQLVKDRLG  
GRSGGYSSGKQGRSDSSSKKDVIELTDDSFDKNVLDSVMMVEFYAPWCGHCCKNLEPEWAAAASEVKEQ  
TKGKVKLAAVDATVNQVLASRYGIRGFPTIKIFQKGESPDYDGGRTSDIVSRALDLFSDNAPPPELLE  
IINEDIKRTCEEHQLCVAVLPHILDTGAAGRNSYLEVLLKLADKYKKMVGWLWTEAGAQSELETALG  
IGGFGYPAMAANARKMKFALLKGSFSEQGINEFLRELSFGRGSTAPVGGGAFPTIVEREPWDGRDGELP  
VEDDIDLSDELDDLGKDEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	47.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_005733</a>



[View online »](#)

Locus ID: 10130

UniProt ID: [Q15084](#), [A0A384NPU5](#)

RefSeq Size: 2349

Cytogenetics: 2p25.1


RefSeq ORF: 1320

Synonyms: ERP5; P5; TXNDC7

**Summary:** This gene encodes a member of the disulfide isomerase (PDI) family of endoplasmic reticulum (ER) proteins that catalyze protein folding and thiol-disulfide interchange reactions. The encoded protein has an N-terminal ER-signal sequence, two catalytically active thioredoxin (TRX) domains, a TRX-like domain, and a C-terminal ER-retention sequence. This protein inhibits the aggregation of misfolded proteins and exhibits both isomerase and chaperone activity. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2016]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified PDIA6 protein (Cat# [TP301710]). The protein was produced from HEK293T cells transfected with PDIA6 cDNA clone (Cat# [RC201710]) using MegaTran 2.0 (Cat# [TT210002]).