

## Product datasheet for **TP301758L**

### DDX23 (NM\_004818) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins  
**Description:** Recombinant protein of human DEAD (Asp-Glu-Ala-Asp) box polypeptide 23 (DDX23), 1 mg  
**Species:** Human  
**Expression Host:** HEK293T  
**Expression cDNA** >RC201758 protein sequence  
**Clone or AA Sequence:** Red=Cloning site Green=Tags(s)

MAGELADKKDRDASPSKEERKRSRTPDRERDRDRDRKSSPSKDRKRHRSRDRRRGGSRSRSRSRSKSAER  
ERRHKERERDKERDRNKKDRDRDKDGHRRDKDRKRSSLSPGRGKDFKSRKDRDSKKDEEDEHGDKKPKAQ  
PLSLEELLAKKKAEEEEAEAKPKFLSKAEREAEALKRRQQEVEERQRMLEEERKKRKQFQDLGRKMLEDPQ  
ERERRERRERMERETNGNEDEEGRQKIREEKDKSKELHAIKERYLGGIKRRRTRHLNDRKFVFEWDASE  
DTSIDYNPLYKERHQVQLLGRGFIAGIDLKQKREQSRFYGDLMEKRRTLEEKEQEEARLRKLRKKEAKQ  
RWDDRHWSSQKKLDEMTDRDWRIFREDYSITTKGGKIPNPIRSWKDSSLPPHILEVIDKCGYKEPTPIQRQ  
AIPIGLQNRDIIGVAETGSGKTA AFLIPLLWITTLPKIDRIEESDQGPYAILAPTRELAQQIEEETIK  
FGKPLGIRTVAVIGGISREDQGFRLRMGCEIATPGRLLIDVLENRYLVLSRCTYVVLDEADRMIDMGFE  
PDVQKILEHMPVSNQKPDTDEAEDPEKMLANFESGKHKYRQTMFTATMPPAVERLARSYLRRPAWYIG  
SAGKPHERVEQKVFLMSESEKRKLLAILEQGFDPPIIFVNQKKGCDVLAKSLEKMGYNACTLHGGKQG  
EQREFALS NLKAGAKDILVATDVAGRGIDIQDVSMVNYDMAKNIEDYIHRIGRTGRAGKSGVAITFLT  
EDSAVFYELKQAILESPVSSCPPPELANHPDAQHKPGTILTKKRREETIFA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Tag:** C-Myc/DDK  
**Predicted MW:** 95.4 kDa  
**Concentration:** >0.05 µg/µL as determined by microplate BCA method  
**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining  
**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol  
**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.  
**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online >](#)

<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_004809</a>
<b>Locus ID:</b>	9416
<b>UniProt ID:</b>	<a href="#">Q9BUQ8</a> , <a href="#">A0A024R0Z3</a> , <a href="#">B3KY11</a>
<b>RefSeq Size:</b>	3288
<b>Cytogenetics:</b>	12q13.12
<b>RefSeq ORF:</b>	2460
<b>Synonyms:</b>	prp28; PRPF28; SNRNP100; U5-100K; U5-100KD
<b>Summary:</b>	This gene encodes a member of the DEAD box protein family. DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure, such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. The protein encoded by this gene is a component of the U5 snRNP complex; it may facilitate conformational changes in the spliceosome during nuclear pre-mRNA splicing. An alternatively spliced transcript variant has been found for this gene, but its biological validity has not been determined. [provided by RefSeq, Jul 2008]
<b>Protein Pathways:</b>	Spliceosome

### Product images:



Coomassie blue staining of purified DDX23 protein (Cat# [TP301758]). The protein was produced from HEK293T cells transfected with DDX23 cDNA clone (Cat# [RC201758]) using MegaTran 2.0 (Cat# [TT210002]).