

Product datasheet for TP301762M

Hsp40 (DNAJB1) (NM_006145) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Dnaj (Hsp40) homolog, subfamily B, member 1 (DNAJB1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201762 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGKDYYQTLGLARGASDEEIKRAYRRQALRYHPDKNKEPGAEEKFKEIAEAYDVLSDPRKREIFDRYGEE GLKGSGPSGGSGGGANGTSFSYTFHGDPHAMFAEFFGGRNPFDTFFGQRNGEEMDIDDPFSGFPMGMGG FTNVNFGRSRQAQEPARKKQDPPVTHDLRVSLEEIYSGCTKKMKISHKRLNPDGKSIRNEDKILTIEVKK GWKEGKITFPKEGDQTSNNIPADIVFLKDKPHNIFKRDGSDVIYPARISLREALCGCTVNVPTLDGRT IPVFKDVIRPGMRRKVPGEGLPLPKTPEKRGDLIEFEVIFPERIPQTSRTVLEQVLP</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	37.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_006136</u>
Locus ID:	3337




[View online »](#)

UniProt ID: [P25685, Q6FHS4](#)
RefSeq Size: 2419
Cytogenetics: 19p13.12
RefSeq ORF: 1020
Synonyms: Hdj1; Hsp40; HSPF1; RSPH16B; Sis1

Summary: This gene encodes a member of the DnaJ or Hsp40 (heat shock protein 40 kD) family of proteins. DNAJ family members are characterized by a highly conserved amino acid stretch called the 'J-domain' and function as one of the two major classes of molecular chaperones involved in a wide range of cellular events, such as protein folding and oligomeric protein complex assembly. The encoded protein is a molecular chaperone that stimulates the ATPase activity of Hsp70 heat-shock proteins in order to promote protein folding and prevent misfolded protein aggregation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]

Product images:



Coomassie blue staining of purified DNAJB1 protein (Cat# [TP301762]). The protein was produced from HEK293T cells transfected with DNAJB1 cDNA clone (Cat# [RC201762]) using MegaTran 2.0 (Cat# [TT210002]).