

Product datasheet for TP301796L

HNRNPD (NM_031370) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human heterogeneous nuclear ribonucleoprotein D (AU-rich element RNA binding protein 1, 37kDa) (HNRNPD), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC201796 representing NM_031370
Clone or AA	Red=Cloning site Green=Tags(s)
Sequence:	<p>MSEEQFGGDGAAAAATAAVGGSAGEQEGAMVAATQGAAAAAGSGAGTGGGTASGGTEGGSAAESEGAKIDA SKNEEDEGHSNSSPRHSEAATAQREEWKFIGGLSWDTTKKDLKDYFSKFGVEVDCTLKLDPITGRSRGF GFVLFKESESVDKVMQKEHKLNGKVIDPKRAKAMKTEPVKKIFVGGGLSPDTPEEKIREYFGGFGEVES IELPMDNKTNRKRGFCFITFKEEEPVKKIMEKKYHNVGLSKCEIKVAMSKEQYQQQQWWSRGGFAGRAR GRGGGPSQWNWQGYSNYWNQGYGNYGYSQGYGGYGGYDYGYNYYGYGDYSNQSGYGVKVSRRGGHQ SYKPY</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	38.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_112738</u>



[View online »](#)

Locus ID: 3184

UniProt ID: [Q14103](#), [A1LU37](#)

RefSeq Size: 2257

Cytogenetics: 4q21.22


RefSeq ORF: 1065

Synonyms: AUF1; AUF1A; hnRNPD0; HNRPD; P37

Summary: This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are nucleic acid binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has two repeats of quasi-RRM domains that bind to RNAs. It localizes to both the nucleus and the cytoplasm. This protein is implicated in the regulation of mRNA stability. Alternative splicing of this gene results in four transcript variants. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified HNRNPD protein (Cat# [TP301796]). The protein was produced from HEK293T cells transfected with HNRNPD cDNA clone (Cat# [RC201796]) using MegaTran 2.0 (Cat# [TT210002]).