

Product datasheet for TP301834M

HNRPH1 (HNRNPH1) (NM_005520) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human heterogeneous nuclear ribonucleoprotein H1 (H) (HNRNPH1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201834 protein sequence Red =Cloning site Green =Tags(s)
	MMLGTEGGEGFVVKVRGLPWSCSADEVQFFSDCKIQNGAQGIRFIYTREGRPSGEAFVELESEDEVKLA LKKDRETMGHRYVEVFKSNNVEMDWLKHGTGPNSPDTANDGFVRLRGLPFGCSKEEIVQFFSGLEIVPN G ITLPVDFQGRSTGEAFVQFASQEIAEKALKKHKERIGHRYIEIFKSSRAEVRTHYDPPRKLMMAMQRPGPY DRPGAGRGYNSIGRGAGFERMRRGAYGGGYGGYDDYNGYNDGYGFGSDRFRDLNYCFSGMSDHRYGD GG STFQSTTGHCVHMRGLPYRATENDIYNFFSPLNPVRVHIEIGPDGRVTGEADVEFATHEDAVAAMSKDKA NMQHRYVELFLNSTAGASGGAYEHRYVELFLNSTAGASGGAYGSQMMGGMGLSNQSSYGGPASQQLS GGY GGGYGGQSSMSGYDQVLQENSSDFQSNIA
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	49 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005511
Locus ID:	3187
UniProt ID:	P31943
RefSeq Size:	2274
Cytogenetics:	5q35.3
RefSeq ORF:	1347
Synonyms:	hnRNPH; HNRPH; HNRPH1
Summary:	<p>This gene encodes a member of a subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA. These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some may shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNA and is very similar to the family member HNRPF. This gene may be associated with hereditary lymphedema type I. Alternatively spliced transcript variants have been described [provided by RefSeq, Mar 2012]</p>

Product images:



Coomassie blue staining of purified HNRNPH1 protein (Cat# [TP301834]). The protein was produced from HEK293T cells transfected with HNRNPH1 cDNA clone (Cat# [RC201834]) using MegaTran 2.0 (Cat# [TT210002]).