

Product datasheet for **TP301870M**

RRNAD1 (NM_015997) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 66 (C1orf66), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201870 protein sequence Red =Cloning site Green =Tags(s)
	MPGISARGLSHEGRKQLAVNLTRVLALYRSILDAYIIFFTDNLWDTLPCSWQEALDGLKPPQLATMLLG MPGEGEVVRYRSVWPLTLLALKSTACALAFTRMPGFQTPSEFLENPSQSSRLTAPFRKHVRPKKQHEIRR LGELVKKLSDFTGCTQVVDVGSGQGHLRFMALGLGLMVKSIEGDQRLVERAQRDQELLQALEKEEKRN PQVWQTSRPHSPHHVVRWVDPTALCEELLPLENPCQGRARLLTGLHACGDLVALLRHFSCCPEVAL ASVGCCYMKLSDPGGYPLSQWVAGLPGYELPYRLREGACHALEEYAERLQKAGPGLRTHCYRAALETVIR RARPELRRPGVQGIPRVHELKIEEYVQRGLQRVGLDPQLPLNLAALQAHLAQENRVVAFFSLALLAPLV ETLILLDRLLYLQEQLFHAELLPIFSPELSPRNLLVLVATKMPLGQALSLETEDS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	52.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_057081</u>
Locus ID:	51093
UniProt ID:	<u>Q96FB5</u>
RefSeq Size:	2283
Cytogenetics:	1q23.1
RefSeq ORF:	1425
Synonyms:	C1orf66; CGI-41; METTL25B

Product images:


Coomassie blue staining of purified RRNAD1 protein (Cat# [TP301870]). The protein was produced from HEK293T cells transfected with RRNAD1 cDNA clone (Cat# [RC201870]) using MegaTran 2.0 (Cat# [TT210002]).