

Product datasheet for **TP301884**

STAMBPL1 (NM_020799) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human STAM binding protein-like 1 (STAMBPL1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201884 protein sequence Red =Cloning site Green =Tags(s)
	<p>MDQPFTVNSLKKLAAMPDHTDVLSPEERVRLSKLGCNITISEDITPRRYFRSGVEMERMASVYLEEGN LENAFVLYNKFITLFVEKLPNHRDYQQCAVPEKQDIMKKLKEIAFPRTDELKNDLLKKYNVEYQEYLQSK NKYKAEILKKLEHQRLIEAERKRIAQMRQQLESEQFLFFEDQLKKQELARGQMRSQQTSGLSEQIDGSA LSCFSTHQNNLLNVFADQPNKSDATNYASHSPVNRALTPAATLSAVQNLVVEGLRCVLPEDLCHKFL QLAESNTVRGIETCGILCGKLTHNEFTITHVIVPKQSAGPDYCDMENVEELFNVQDQHDLLTLGWIHTHP TQTAFLSSVDLHHTCSYQLMLPEAIAIVCSPKHKDTGIFRLTNAGMLEVSACKKKGFHPHTKEPRLFSIC KHVLVKDIKIIVLCLR</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	49.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_065850</u>



[View online »](#)

Locus ID:	57559
UniProt ID:	Q96FJ0
RefSeq Size:	2094
Cytogenetics:	10q23.31
RefSeq ORF:	1308
Synonyms:	ALMalpha; AMSH-FP; AMSH-LP; bA399O19.2
Summary:	Zinc metalloprotease that specifically cleaves 'Lys-63'-linked polyubiquitin chains. Does not cleave 'Lys-48'-linked polyubiquitin chains.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Protease

Product images:



Coomassie blue staining of purified STAMBPL1 protein (Cat# TP301884). The protein was produced from HEK293T cells transfected with STAMBPL1 cDNA clone (Cat# [RC201884]) using MegaTran 2.0 (Cat# [TT210002]).