

## Product datasheet for TP301902M

### FAM118A (NM\_017911) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 118, member A (FAM118A), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201902 protein sequence Red=Cloning site Green=Tags(s)
	<p>MDSVEKTTNRSEQKSRKFLKSLIRKQPQELLLVIGTGVSAAVAPGIPALCSWRSCIEAVIEAAEQLEVLH PGDVAEFRKVKTKDRDLLVVAHDLIRKMSPRTGDAKPSFFQDCLMEVFDDLEQHRSPLVLQSILSLMDR GAMVLTTNYDNLLEAFGRQNKPMESLDLKDQTKVLEWARGHMKYGVLHIHGLYRDPGCVLDPSPGYKDV TQDAEVMEVLQNLRYRKSFLVGCGETLHDQIFQALFLYVSNKVDLEHYMLVLKENEDHFFKHQADMLL HGKVVSYGDCFDHFPQVQDLATQICKQQSPDADRVDSTLLGNACQDCAKRKLENGIEVSKKRTQSD TDDAGGS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	40.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_060381</a></u>



[View online »](#)

Locus ID: 55007  
UniProt ID: [Q9NWS6](#)  
RefSeq Size: 2928  
Cytogenetics: 22q13.31  
RefSeq ORF: 1071  
Synonyms: C22orf8

**Product images:**

Coomassie blue staining of purified FAM118A protein (Cat# [TP301902]). The protein was produced from HEK293T cells transfected with FAM118A cDNA clone (Cat# [RC201902]) using MegaTran 2.0 (Cat# [TT210002]).