

## Product datasheet for TP301925L

### 14 3 3 gamma (YWHAG) (NM\_012479) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, gamma polypeptide (YWHAG), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201925 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MVDREQLVQKARLAEQAERYDDMAAAMKNVTELNEPLSNEERNLLSVAYKNVVGARRSSWRVISSIEQKT SADGNEKKIEMVRAYREKIEKELEAVCQDVLSDNYLIKNCSETQYESKVFYLMKMGDYRYLAEVATG EKRA TVESSEKAYSEAHEISKEHMQPTHPIRLGLALNYSVFYIEIQNAPEQACHLAKTAFDDAIAELDT LNESYKDSLIMQLLRDNLTLWTSQQDDGGEGNN  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span>
Tag:	C-Myc/DDK
Predicted MW:	28.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_036611</a></u>
Locus ID:	7532


[View online »](#)

UniProt ID: [P61981](#)

RefSeq Size: 3779

Cytogenetics: 7q11.23

RefSeq ORF: 741


Synonyms: 14-3-3GAMMA; DEE56; EIEE56; PPP1R170

**Summary:** This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 100% identical to the rat ortholog. It is induced by growth factors in human vascular smooth muscle cells, and is also highly expressed in skeletal and heart muscles, suggesting an important role for this protein in muscle tissue. It has been shown to interact with RAF1 and protein kinase C, proteins involved in various signal transduction pathways. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis

### Product images:



Coomassie blue staining of purified YWHAG protein (Cat# [TP301925]). The protein was produced from HEK293T cells transfected with YWHAG cDNA clone (Cat# [RC201925]) using MegaTran 2.0 (Cat# [TT210002]).