

Product datasheet for **TP301930**

TRIM31 (NM_007028) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tripartite motif-containing 31 (TRIM31), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201930 representing NM_007028 Red =Cloning site Green =Tags(s) MASGQFVNKLQEEVICPICLDILQKPVTIDCGHNFCLKCITQIGETSCGFFKCPLCKTSVRKNAIRFNSL LRNLVEKIQALQASEVQSKRKEATCPRHQEMFHYFCEDDGKFLCFVCRESKDHKSHNVSLIEEAAQNYQG QIQEQIQLQQKEKETVQVKAQGVHRVDVFTDQVEHEKQRILTEFELLHQVLEEEKNFLLSRIYWLGHG TEAGKHYVASTEPQLNDLKKLVDSLKTKQNMPPRQLLEDIKVVLCSRSEEFQLNPTVPVPLEKKLSEAK SRHDSITGSLKKFKDQLQADRKKDENRFFKSMNKNDMKSWGLLQKNNHKMNKTSEPGSSSAGGRTTSG PP NHHSSAPSHSLFRASSAGKVTFPVCLLASYSDEISGQGASSQDTKTFDVALSEELHAALSEWLTAIRAWFC KVPSS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	48.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_008959</u>
Locus ID:	11074
UniProt ID:	<u>Q9BZY9</u>
RefSeq Size:	2049
Cytogenetics:	6p22.1
RefSeq ORF:	1275
Synonyms:	C6orf13; HCG1; HCGI; RNF
Summary:	This gene encodes a protein that functions as an E3 ubiquitin-protein ligase. This gene shows altered expression in certain tumors and may be a negative regulator of cell growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]

Product images:



Coomassie blue staining of purified TRIM31 protein (Cat# TP301930). The protein was produced from HEK293T cells transfected with TRIM31 cDNA clone (Cat# [RC201930]) using MegaTran 2.0 (Cat# [TT210002]).