

Product datasheet for TP301994L

XPG (ERCC5) (NM_000123) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human excision repair cross-complementing rodent repair deficiency, complementation group 5 (ERCC5), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201994 protein sequence Red=Cloning site Green=Tags(s)

MGVQGLWKLLCESGRQVSPEALEGKILAVDISIWLNQALKGVRDRHGNSIENPHLLTLFHRLCKLLFFRI
RPIFVFDGDAPLLKKQTLVKRRQRKDLASSDSRKTTEKLLKTFKLRQAIKTAFRSKRDEALPSLTQVRRE
NDLYVLPPLQEEKHSSSEEEDEKEWQERMNQKALQEEFFHNPQAIDIESEDFSSLPPEVKHEILTDMKE
FTKRRRTLFEAMPEESDDFSQYQLKGLLKKNYLNQHIEHVQKEMNQQHSGHIRRQYEDEGGFLKEVESRR
VSEDTSHYILIKGIQAKTVAEVDSESLPSSSKMHGMSFDVKSSPCEKLTEKEPDATPPSPRTLLAMQA
ALLGSSSEEELESENRRQARGRNAPAAVDEGSISPRTLAIKRALDDDEDVKVCAGDDVQTGGPGAEMR
INSSTENSDEGLKVRDGGKIPFTATLASSSVNSAEHVASTNEGREPTDVPKEQMSLVHVGTEAFPISD
ESMIKDRKDRPLPESAVVRHSDAPGLPNGRELTPASPTCTNSVSKNETHAEVLEQQNELCPYESKFDSSL
LSSDDETKCKPNSASEVIGPVSLQETSSIVSVPSEAVDNVENVVSFNAKEHENFLETIQEQTTESAGQD
LISIPKAVEPMEIDSEESDGSFIEVQSVISDEELQAEFPETSKPPSEQGEEELVGTREGEAPAESESL
LRDNSERDDVDGEPQAEKDAEDSLHEWQDINLEELETLESNLLAQQNSLKAQKQQQERIAATVTGQMFL
ESQELLRLFGIPYIQAPMEAEAQCAILDLTDTQSGTITDDSDIWLFGARHVYRNFFNKNKFVEYYQYVDF
HNQLGLDRNKLINLAYLLGSDYTEGIPTVGCVTAMEILNEFPGHGLEPLKLFSEWWHEAQKNPKIRPNPH
DTKVKKKLRTLQLTPGFNPVAEAYLKPVVDDSKGSFLWGKPDLDKIREFCQRYFGWNRTKTDESFPV
LKQLDAQQTQLRIDSFFRLAQQEKEDAKRIKSQRLNRAVTCMLRKEKEAAAASEIAVSVAMEKEFELLDK
AKRKTQKRGITNTLESSSLKRKRLSDSKRKNCTCGGFLGETCLSESSDGSSEDAEASSSLMNVQRRTAAK
EPKTSASDSQNSVKEAPVKNGGATTSSSSDSDDDGGKEMVLVTARSVFGKRRRKLRRARGRKRKT


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	133.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining



[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000114
Locus ID:	2073
UniProt ID:	P28715
RefSeq Size:	4091
Cytogenetics:	13q33.1
RefSeq ORF:	3558
Synonyms:	COFS3; ERCC5-201; ERCM2; UVDR; XPG; XPGC
Summary:	This gene encodes a single-strand specific DNA endonuclease that makes the 3' incision in DNA excision repair following UV-induced damage. The protein may also function in other cellular processes, including RNA polymerase II transcription, and transcription-coupled DNA repair. Mutations in this gene cause xeroderma pigmentosum complementation group G (XP-G), which is also referred to as xeroderma pigmentosum VII (XP7), a skin disorder characterized by hypersensitivity to UV light and increased susceptibility for skin cancer development following UV exposure. Some patients also develop Cockayne syndrome, which is characterized by severe growth defects, cognitive disability, and cachexia. Read-through transcription exists between this gene and the neighboring upstream BIVM (basic, immunoglobulin-like variable motif containing) gene. [provided by RefSeq, Feb 2011]
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transcription Factors
Protein Pathways:	Nucleotide excision repair

Product images:

Coomassie blue staining of purified ERCC5 protein (Cat# [TP301994]). The protein was produced from HEK293T cells transfected with ERCC5 cDNA clone (Cat# [RC201994]) using MegaTran 2.0 (Cat# [TT210002]).