

Product datasheet for **TP302000**

CD9 (NM_001769) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human CD9 molecule (CD9), 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC202000 protein sequence
Red=Cloning site **Green**=Tags(s)

MPVKGGTKCIKYLFLGFNFIFWLAGIAVLAIGLWLRFDSQTKSIFEQETNNNNSSFYTGVIILIGAGALM
 MLVGFLGCCGAVQESQCMLGLFFGFLLVIFAIEIAAAIWGYSHKDEVIKEVQEFYKDTYNKLTDEPQR
 ETLKAIHYALNCCGLAGGVEQFISDICPKKDVLETFTVKSCPDAIKEVFDNKFHIIIGAVGIGIAVVMIFG
 MIFSTILCCAIRRNREMV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 25.2 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity: ELISA positive control (PMID: [26857499](https://pubmed.ncbi.nlm.nih.gov/26857499/))
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_001760](https://www.ncbi.nlm.nih.gov/RefSeq/record/NP_001760)
Locus ID: 928



[View online »](#)

UniProt ID: [P21926](#)

RefSeq Size: 1321

Cytogenetics: 12p13.31

RefSeq ORF: 684


Synonyms: BTCC-1; DRAP-27; MIC3; MRP-1; TSPAN-29; TSPAN29

Summary: This gene encodes a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Tetraspanins are cell surface glycoproteins with four transmembrane domains that form multimeric complexes with other cell surface proteins. The encoded protein functions in many cellular processes including differentiation, adhesion, and signal transduction, and expression of this gene plays a critical role in the suppression of cancer cell motility and metastasis. [provided by RefSeq, Jan 2011]

Protein Families: Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transmembrane

Protein Pathways: Hematopoietic cell lineage

Product images:



Coomassie blue staining of purified CD9 protein (Cat# TP302000). The protein was produced from HEK293T cells transfected with CD9 cDNA clone (Cat# [RC202000]) using MegaTran 2.0 (Cat# [TT210002]).