

Product datasheet for TP302026

SCO1 (NM_004589) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SCO cytochrome oxidase deficient homolog 1 (yeast) (SCO1), nuclear gene encoding mitochondrial protein, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202026 protein sequence Red=Cloning site Green=Tags(s)

MAMLVLVPGRVMRPLGGQLWRFLPRGLEFWGPAEGTARVLLRQFCARQAEAWRASGRPGYCLGTRPLSTA
RPPPPWSQKGGDSTRPSKPGPVSWKSLAITFAIGGALLAGMKHVKKEKAEKLEKERQRHIGKPLLGGPF
SLTHTGERKTDKDYLGQWLLIYFGFTHCPDVCPEELEKMIQVVDEIDSITLPLDLPLFISIDPERDTK
EAIANYVKEFSPKLVGLTGTREEVDQVARAYRVVYSPGPKDEDEDYIVDHTIIMYLIGPDGFLDYFGQN
KRGGEIAASIATHMRPYRKKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	33.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004580
Locus ID:	6341



[View online »](#)

UniProt ID:	O75880
RefSeq Size:	1768
Cytogenetics:	17p13.1
RefSeq ORF:	903
Synonyms:	MC4DN4; SCOD1
Summary:	Mammalian cytochrome c oxidase (COX) catalyzes the transfer of reducing equivalents from cytochrome c to molecular oxygen and pumps protons across the inner mitochondrial membrane. In yeast, 2 related COX assembly genes, SCO1 and SCO2 (synthesis of cytochrome c oxidase), enable subunits 1 and 2 to be incorporated into the holoprotein. This gene is the human homolog to the yeast SCO1 gene. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified SCO1 protein (Cat# TP302026). The protein was produced from HEK293T cells transfected with SCO1 cDNA clone (Cat# [RC202026]) using MegaTran 2.0 (Cat# [TT210002]).