

Product datasheet for **TP302071M**

Noxa (PMAIP1) (NM_021127) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phorbol-12-myristate-13-acetate-induced protein 1 (PMAIP1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202071 representing NM_021127 Red =Cloning site Green =Tags(s)
	MPGKKARKNAQPSPARAPAELEVECATQLRRFGDKLNFRQKLLNLISKLFCSGT
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	5.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_066950
Locus ID:	5366
UniProt ID:	Q13794
RefSeq Size:	1885
Cytogenetics:	18q21.32



[View online »](#)

RefSeq ORF: 162


Synonyms: APR; NOXA

Summary: This gene belongs to a pro-apoptotic subfamily within the BCL-2 protein family, referred to as the BCL-2 homology domain 3 (BH3)-only subfamily, which determine whether a cell commits to apoptosis. In response to death-inducing stimuli, BH3-only members inhibit the anti-apoptotic BCL-2 family members, which under steady-state conditions keep the multi-BH domain proteins BAX and BAK, in an inactive state. [provided by RefSeq, Aug 2020]

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: p53 signaling pathway

Product images:



Coomassie blue staining of purified PMAIP1 protein (Cat# [TP302071]). The protein was produced from HEK293T cells transfected with PMAIP1 cDNA clone (Cat# [RC202071]) using MegaTran 2.0 (Cat# [TT210002]).