

Product datasheet for TP302093

FA2H (NM_024306) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fatty acid 2-hydroxylase (FA2H), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202093 protein sequence Red =Cloning site Green =Tags(s)
	MAPAPPPAASFSPSEVQRRRLAAGACWVRRGARLYDLSSFVRHHPGGEQLLRARAGQDISADLDGPPHRHS ANARRWLEQYYVGELRGEQQGSMENEPVALEETQKTDPAPEPRFKVVDWDKDLVDWRKPLLWQVGHGK YDEWVHQPVTRPIRLFHSDLIEGLSKTVWYSVPIIWVPLVLYLSWSYYRTFAQGNVRLFTSFTTEYTVAV PKSMFPGFLMGTFLWSLIEYLIHRFLFHMKPPSDSYLIMLHFVMHGQHHKAPFDGSRVFPVPASLV IGVFYLCMQLILPEAVGGTVFAGGLLGYVLYDMTHYYLHFGSPHKGSYLYSLKAHHVKHHFAHQSGFGI STKLWDYCFHTLTPEKPHLKTQ
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	42.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_077282
Locus ID:	79152




[View online »](#)

UniProt ID: [Q7L5A8](#)
RefSeq Size: 2450
Cytogenetics: 16q23.1
RefSeq ORF: 1116
Synonyms: FAAH; FAH1; FAXDC1; SCS7; SPG35

Summary: This gene encodes a protein that catalyzes the synthesis of 2-hydroxysphingolipids, a subset of sphingolipids that contain 2-hydroxy fatty acids. Sphingolipids play roles in many cellular processes and their structural diversity arises from modification of the hydrophobic ceramide moiety, such as by 2-hydroxylation of the N-acyl chain, and the existence of many different head groups. Mutations in this gene have been associated with leukodystrophy dysmyelinating with spastic paraparesis with or without dystonia.[provided by RefSeq, Mar 2010]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified FA2H protein (Cat# TP302093). The protein was produced from HEK293T cells transfected with FA2H cDNA clone (Cat# [RC202093]) using MegaTran 2.0 (Cat# [TT210002]).