

Product datasheet for **TP302130M**

GSTA4 (NM_001512) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human glutathione S-transferase alpha 4 (GSTA4), 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC202130 protein sequence
Red=Cloning site **Green**=Tags(s)

MAARPKLHYPNGRGRMESVRWVLAAGVEFDEEFLETKEQLYKLQDGNHLLFQQVPMVEIDGMKLVQTRS
ILHYIADKHNLFGKNLKERTLIDMYVEGTLDLLELLIMHPFLKPDDQQKEVWNMAQKAIIRYFPVFEKIL
RGHGQSFLVGNQLSLADVILLQTILALEEKIPNILSAFPFLQEYTVKLSNIPTIKRFLEPGSKKKPPPDE
IYVRTVYNIFRP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 25.5 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001503](#)

Locus ID: 2941

UniProt ID: [O15217](#), [A0A024RD58](#)



[View online >](#)

RefSeq Size: 1352

Cytogenetics: 6p12.2


RefSeq ORF: 666

Synonyms: GSTA4-4; GTA4

Summary: Cytosolic and membrane-bound forms of glutathione S-transferase are encoded by two distinct supergene families. These enzymes are involved in cellular defense against toxic, carcinogenic, and pharmacologically active electrophilic compounds. At present, eight distinct classes of the soluble cytoplasmic mammalian glutathione S-transferases have been identified: alpha, kappa, mu, omega, pi, sigma, theta and zeta. This gene encodes a glutathione S-transferase belonging to the alpha class. The alpha class genes, which are located in a cluster on chromosome 6, are highly related and encode enzymes with glutathione peroxidase activity that function in the detoxification of lipid peroxidation products. Reactive electrophiles produced by oxidative metabolism have been linked to a number of degenerative diseases including Parkinson's disease, Alzheimer's disease, cataract formation, and atherosclerosis. [provided by RefSeq, Jul 2008]

Protein Pathways: Drug metabolism - cytochrome P450, Glutathione metabolism, Metabolism of xenobiotics by cytochrome P450

Product images:



Coomassie blue staining of purified GSTA4 protein (Cat# [TP302130]). The protein was produced from HEK293T cells transfected with GSTA4 cDNA clone (Cat# [RC202130]) using MegaTran 2.0 (Cat# [TT210002]).