

## Product datasheet for TP302162

### LIMPII (SCARB2) (NM\_005506) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human scavenger receptor class B, member 2 (SCARB2), 20 µg

**Species:** Human

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >RC202162 protein sequence  
Red=Cloning site Green=Tags(s)

MGRCCFYTAGTLSLLLLVTSVTLVARVFQKAVDQSIKIVLRNGTEAFDSWEKPPLPVYTQFYFFNVT  
NPEEILRGETPRVEEVGPYTYRELRNKANIQFGDNGTTISAVSNKAYVFERDQSVGDPKIDLIRTLNIPV  
LTVIEWSQVHFLREIIEAMLKAYQQKLFVTHTVDELLWGYKDEILSLIHVFRPDISPYGLFYEKNGTND  
GDYVFLTGEDSYLNFTKIVEWNGKTSLDWWITDKCNMINGTDGDSFHPLITKDEVLYVFPDFCRSVYIT  
FSDYESVQGLPAFRYKVP AEILANTSDNAGFCIPEGNCLGSGVLNVSICKNGAPIIMSFPHFYQADERFV  
SAIEGMHPNQEDHETFVDINPLTGILKAAKRFQINIYVKKLDDFVETGDIRTMVFPV MYLNESVHIDKE  
TASRLKSMINTTLIITNIPYIIMALGVFFGLVFTWLACKGQGSMDDEGTADERAPLIRT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Tag:** C-Myc/DDK

**Predicted MW:** 54.1 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_005497](#)



[View online »](#)


Locus ID:	950
UniProt ID:	<a href="#">Q14108</a> , <a href="#">A0A024RDG6</a>
RefSeq Size:	4780
Cytogenetics:	4q21.1
RefSeq ORF:	1434
Synonyms:	AMRF; CD36L2; EPM4; HLG85; LGP85; LIMP-2; LIMPII; SR-BII

**Summary:** The protein encoded by this gene is a type III glycoprotein that is located primarily in limiting membranes of lysosomes and endosomes. Earlier studies in mice and rat suggested that this protein may participate in membrane transportation and the reorganization of endosomal/lysosomal compartment. The protein deficiency in mice was reported to impair cell membrane transport processes and cause pelvic junction obstruction, deafness, and peripheral neuropathy. Further studies in human showed that this protein is a ubiquitously expressed protein and that it is involved in the pathogenesis of HFMD (hand, foot, and mouth disease) caused by enterovirus-71 and possibly by coxsackievirus A16. Mutations in this gene caused an autosomal recessive progressive myoclonic epilepsy-4 (EPM4), also known as action myoclonus-renal failure syndrome (AMRF). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2011]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Lysosome

### Product images:



Coomassie blue staining of purified SCARB2 protein (Cat# TP302162). The protein was produced from HEK293T cells transfected with SCARB2 cDNA clone (Cat# [RC202162]) using MegaTran 2.0 (Cat# [TT210002]).