

Product datasheet for TP302168

Calpastatin (CAST) (NM_001042444) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calpastatin (CAST), transcript variant 10, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202168 protein sequence Red =Cloning site Green =Tags(s)

MNPTETKAVKTEPEKKSQSTKPKSLPKQASDTGSNDAHNNKAVSRSAEQQPSEKSTEPKTKPQDMISAGG
ESVAGITAISGKPGDKKKEKSLTPAVPVESKPKPSGKSGMDAALDDLIDTLGGPEETEENTTYTGPE
VSDPMSSTYIEELGKREVTIPPKYRELLAKKEGITGPPADSSKPIGPDDAIDALSSDFTCGSPTAAGKKT
EKEESTEVLKAQSAGTVRSAAPPQEKRRKVEKDTMSDQALEALSASLGTRQAEPELDRSIKEVDEAKAK
EEKLEKCGEDDETIPSEYRLKPATDKDGKPLLPEPEEKPKPRSESELIDELSEDFRSECKEKPSKPTK
TEESKAAAPAPVSEAVCRTSMCSIQSAPPEPATLKGTVPDDAVEALADSLGKKEADPEDGKPVMDKVKEK
AKEEDREKLGEKEETIPPDYRLEEVKDKDGKPLLPKESKEQLPPMSDFLLDALSEDFSGPQNASSLKFE
DAKLA AISEVVSQTPASTTQAGAPPRDTSQSDKLDLDDALDKLSDSLGRQPDNDENKPMEDKVKEKAKA
EHRDKLGERDDTIPPEYRHLLDDNGQDKPVKPPTKKSSEDSKPPADDQDPIDALSGDLSDSCPSTTETSQNT
AKDKCKKAASSSKAPKNGGKAKDSAKTTEETSKPKDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	71.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001035909
Locus ID:	831
UniProt ID:	P20810
RefSeq Size:	4366
Cytogenetics:	5q15
RefSeq ORF:	2001
Synonyms:	BS-17; PLACK
Summary:	<p>The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]</p>

Product images:



Coomassie blue staining of purified CAST protein (Cat# TP302168). The protein was produced from HEK293T cells transfected with CAST cDNA clone (Cat# [RC202168]) using MegaTran 2.0 (Cat# [TT210002]).