

## Product datasheet for TP302186L

### SH3KBP1 (NM\_031892) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SH3-domain kinase binding protein 1 (SH3KBP1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202186 representing NM_031892 Red=Cloning site Green=Tags(s)

MVEAIVEFDYQAQHDELTISVGEIITNIRKEDGGWWEGQINGRRGLFPDNFVREIKKEMKKDPLTNKAP  
EKPLHEVPSGNSLLSSETILRTNKRGERRRRCQVAFSYLPQNDDELELKVGDIIIEVVGEVEEGWWEGVL  
NGKTGMFSPNFIKELSGESDELGISQDEQLSKSSLRETTGSESDGGDSSTKSEGANVTATAIQPKKV  
KGVGFGDIFKDKPIKLRPRSIEVENDFLPVEKTIGKKLPATTATPDSSKTEMDSRTKSKDYCKVIFPYEA  
QNDDELTIKEGDIVTLINKDCIDVGWWEGELNGRRGVFPDNFVKLLPPDFEKEGNRPKPPPPSAPVIKQ  
GAGTTERKHEIKKIPPERPEMLPNRTEEKERPEREKLDLQKPSVPAIPPKPRPPKTNLSLRPGALPPR  
RPERPVGPLTHTRGDSPKIDLAGSSLSGILDKDLSDRSNDIDLEGFDSVVSSTEKLSHPTTSRPKATGRR  
PPSQSLTSSSLSPDIFDPSPEEDKEEHISLAHRGVDASKKTSKTVTISQVSDNKASLPPKPGTMAAGG  
GGPAPLSSAVPSPLSSSLGTAGHRANSPSLFGTEGKPKMEPAASSQAAVEELRTQVRELRSIETMKDQQ  
KREIKQLSELDEEKIRLRLQMEVNDIKKALQSK

TRRLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	72.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_114098</a>
<b>Locus ID:</b>	30011
<b>UniProt ID:</b>	<a href="#">Q96B97</a> , <a href="#">Q5JPT6</a>
<b>RefSeq Size:</b>	3348
<b>Cytogenetics:</b>	Xp22.12
<b>RefSeq ORF:</b>	1995
<b>Synonyms:</b>	AGMX2; CD2BP3; CIN85; GIG10; HSB-1; HSB1; IMD61; MIG18
<b>Summary:</b>	This gene encodes an adapter protein that contains one or more N-terminal Src homology domains, a proline rich region and a C-terminal coiled-coil domain. The encoded protein facilitates protein-protein interactions and has been implicated in numerous cellular processes including apoptosis, cytoskeletal rearrangement, cell adhesion and in the regulation of clathrin-dependent endocytosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2017]
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Endocytosis

### Product images:



Coomassie blue staining of purified SH3KBP1 protein (Cat# [TP302186]). The protein was produced from HEK293T cells transfected with SH3KBP1 cDNA clone (Cat# [RC202186]) using MegaTran 2.0 (Cat# [TT210002]).