

Product datasheet for TP302201

VAT1 (NM_006373) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human vesicle amine transport protein 1 homolog (T. californica) (VAT1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202201 protein sequence Red=Cloning site Green=Tags(s)

MSDEREVAEAATGEDASSPPPKTEAASDPQHHPAASEGAAAAAASPPLLRLCLVLTGFGGYDKVKLQSRPAA
PPAPGPGQLTLRLRACGLNFADLMARQGLYDRLPPLPVTTPGMEGAGVVIAGGEGVSDRKAGDRVMVLNRS
GMWQEEVTVPSVQTFLLIPEAMTFEEAAALLVNYITAYMVLDFGNLQPGHSLVHMAAGGVMMAAVQLCR
TVENVTVFGTASASKHEALKENGVTHTPIDYHTTDYVDEIKKISPKGVDIVMDPLGGSDTAKGYNLLKPMG
KVTYGMANLLTGPKRNLMALARTWWNQFSVTALQLLQANRAVCGFHLGYLDGEVELVSGVVARLLALYN
QGHKPHIDSVWPFQKADAMKQMQEKKNVGKVVLLVPGPEKEN

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	41.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006364



[View online »](#)

Locus ID: 10493

UniProt ID: [Q99536](#), [A0A024R1Z6](#)

RefSeq Size: 2758

Cytogenetics: 17q21.31


RefSeq ORF: 1179

Synonyms: VAT1

Summary: Synaptic vesicles are responsible for regulating the storage and release of neurotransmitters in the nerve terminal. The protein encoded by this gene is an abundant integral membrane protein of cholinergic synaptic vesicles and is thought to be involved in vesicular transport. It belongs to the quinone oxidoreductase subfamily of zinc-containing alcohol dehydrogenase proteins. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified VAT1 protein (Cat# TP302201). The protein was produced from HEK293T cells transfected with VAT1 cDNA clone (Cat# [RC202201]) using MegaTran 2.0 (Cat# [TT210002]).