

Product datasheet for **TP302209L**

Hsc70 (HSPA8) (NM_006597) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human heat shock 70kDa protein 8 (HSPA8), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202209 protein sequence Red =Cloning site Green =Tags(s)

MSKGPVAVGIDLGTTYSCVGVFQHGKVEIIANDQGNRTTPSYVAFTDTERLIGDAAKNQVAMNPTNTVFDA
KRLIGRRFDDAVVQSDMKHWPFMVNDAGRPKVQVEYKGETKSFYPEEVSSMVLTKMKEIAEAYLGKTVT
NAVVTVPAYFNDSQRQATKDAGTIAGLNVLRIINEPTAAAIAYGLDKKVGAEARNVLIFFDLGGGTFDVSIL
TIEDGIFEVKSTAGDTHLGGEDFDNRMVNHFAIEFKRKHKKDISENKRAVRLRTACERAKRTLSSSTQA
SIEIDSLYEGIDFYTSITRARFEELNADLFRGTLDPVEKALRDAKLDKSQIHDIIVLVGGSTRIPKIQKLL
QDFFNKELNKSINPDEAVAYGAAVQAAILSGDKSENVQDLLLLDVTPLSLGIETAGGVMTVLIKRNTTI
PTKQTQFTTYSNDNQPGLIQVYEGERAMTKDNNLLGKFELTGIPPAPRGVPQIEVTFDIDANGILNVSA
VDKSTGKENKITITNDKGRLSKEDIERMVQEAKEYKADEKQRDKVSSKNSLESYAFNMKATVEDEKLQG
KINDEDKQKILDKCNEIINWLDKNQTAEEFEHQQKELEKVCNPIITKLYQSAGGMPGGMPGGFPGGGA
PPSGGASSGPTIEEVD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	70.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006588
Locus ID:	3312
UniProt ID:	P11142 , V9HW22 , Q53HF2
RefSeq Size:	2473
Cytogenetics:	11q24.1
RefSeq ORF:	1938
Synonyms:	HEL-33; HEL-S-72p; HSC54; HSC70; HSC71; HSP71; HSP73; HSPA10; LAP-1; LAP1; NIP71
Summary:	This gene encodes a member of the heat shock protein 70 family, which contains both heat-inducible and constitutively expressed members. This protein belongs to the latter group, which are also referred to as heat-shock cognate proteins. It functions as a chaperone, and binds to nascent polypeptides to facilitate correct folding. It also functions as an ATPase in the disassembly of clathrin-coated vesicles during transport of membrane components through the cell. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	Antigen processing and presentation, Endocytosis, MAPK signaling pathway, Spliceosome

Product images:

Coomassie blue staining of purified HSPA8 protein (Cat# [TP302209]). The protein was produced from HEK293T cells transfected with HSPA8 cDNA clone (Cat# [RC202209]) using MegaTran 2.0 (Cat# [TT210002]).