

Product datasheet for **TP302215L**

ZMAT5 (NM_019103) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger, matrin type 5 (ZMAT5), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC202215 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MGKRYFCDYCDRSFQDNLHNRKKHLNGLQHLKAKKWWYDMFRDAAAILLDEQNKRPCRKFLLTGQCDFGS
NCRFSHMSEKDLQELSIQVEEERRAREWLLDAPELPEGHLEDWLEKRAKRLSSAPSSRAEPIRTTVFYYP
VGWPPVQELPPSLRAPPDGGWPLQPRVQWG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	19.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_061976
Locus ID:	55954
UniProt ID:	Q9UDW3 , A0A024R111
RefSeq Size:	1051



[View online »](#)

Cytogenetics: 22q12.2
RefSeq ORF: 510
Synonyms: SNRNP20; U11/U12-20K; ZC3H19

Product images:

Coomassie blue staining of purified ZMAT5 protein (Cat# [TP302215]). The protein was produced from HEK293T cells transfected with ZMAT5 cDNA clone (Cat# [RC202215]) using MegaTran 2.0 (Cat# [TT210002]).