

Product datasheet for **TP302229M**

HCE (RNGTT) (NM_003800) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RNA guanylyltransferase and 5'-phosphatase (RNGTT), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202229 representing NM_003800 Red =Cloning site Green =Tags(s)

MAHNKIPPRWLNCPRRGQPVAGRFLPLKTM LGPRYDSQVAEENRFHPSMLSNYLKSLKVKMGLLVDLTNT
SRFYDRNDIEKEGIKYIKLQCKGHGECPTTENTETFIRLCERFNERNPPELIGVHCTHGFNRTGFLICAF
LVEKMDWSIEAAVATFAQARPPGIYKGDYKELFRRYGDIEEAPPPPLLPDWC FEDDEDEDEDEDGKKES
EPGSSASFGKRRKERLKLGAIFLEGVTVKGVTVTTQPKLGEVQKCHQFCGWEGSGFPGAQPVSMQKQ
IKLLDLKPYKVSWKADGTRYMMLIDGTNEVF MIDRDNVSVFVSNLEFPFRKDLRMHLSNTLLDGEMIIDR
VNGQAVPRYLIYDIIFNSQPVGDCDFNVRLQCIEREIISPRHEKMKTGLIDKTQEPFSVRNKPFFDICT
SRKLLGNFAKEVSHEMDGLIFQPTGKYKPGRCDDILKWKPPSLNSVDFRLKITRMGGEGLLPQNVGLLY
VGGYERPFQIKVTKELKQYDNKIIECKFENNSWVFMQRRTDKSFPNAYNTAMAVCNISINPVTKEMLFE
FIDRCTAASQGQKRKHHLDPDELMPPPPPKRPHPLT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	68.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003791
Locus ID:	8732
UniProt ID:	O60942 , Q7Z3R6
RefSeq Size:	4460
Cytogenetics:	6q15
RefSeq ORF:	1791
Synonyms:	CAP1A; hCAP; HCE; HCE1
Summary:	Bifunctional mRNA-capping enzyme exhibiting RNA 5'-triphosphatase activity in the N-terminal part and mRNA guanylyltransferase activity in the C-terminal part. Catalyzes the first two steps of cap formation: by removing the gamma-phosphate from the 5'-triphosphate end of nascent mRNA to yield a diphosphate end, and by transferring the gmp moiety of GTP to the 5'-diphosphate terminus.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Phosphatase

Product images:



Coomassie blue staining of purified RNGTT protein (Cat# [TP302229]). The protein was produced from HEK293T cells transfected with RNGTT cDNA clone (Cat# [RC202229]) using MegaTran 2.0 (Cat# [TT210002]).