

Name: Recombinant protein of human polymerase (RNA) I polypeptide D, 16kDa (POLR1D), transcript variant 2

Catalog: TP302300

Product Data Sheet

Gene Name: Homo sapiens polymerase (RNA) I polypeptide D, 16kDa (POLR1D), transcript variant 2

GenBank accession: NM_152705

Description: Recombinant protein of human polymerase (RNA) I polypeptide D, 16kDa (POLR1D), transcript variant 2

Protein Gel Analysis Data:

Amount: 20 ug

Buffer and Storage: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.

Source: Recombinant protein was produced with TrueORF clone, [RC202300](#), encoding the full-length human POLR1D with C-terminal DDK tag, from human HEK293 cells.

Protein Accn: NP_689918


Gene Synonym: AC19; POLR1C; RPA16; RPA9; RPAC2; RPC16; RPO1-3; TCS2

Predicted molecular weight: 14.2 kDa

Tags: C-terminal MYC/DDK

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining.

Concentration: >50 ug/mL as determined by microplate BCA method



Purified recombinant protein POLR1D was analyzed by SDS-PAGE gel and Coomassie Blue Staining.