

Product datasheet for **TP302314**

BRF2 (NM_018310) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human BRF2, subunit of RNA polymerase III transcription initiation factor, BRF1-like (BRF2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202314 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPGRGRCPDCGSTELVEDSHYSQSQLVCSDCGCVVTEGVLTTFSDDEGNLREVTYSRSTGENEQVSRSQQ RGLRRVRDLCRVLQLPPTFEDTAVAYYQQAYRHSGIRAARLQKKEVLVGCCVLITCRQHNWPLTMGAICT LLYADLDVFSSTYMQIVKLLGLDVPSSLCLAEVLKTYCSSFKLFQASPSVPAKYVEDKEKMLSRMTMQLVEL ANETWLVTGRHPLPVITAATFLAWQSLQPADRLSCSLARFCKLANVDLPYPASSRLQELLAVLLRMAEQL AWLRVLRDLKRSVVKHIGDLLQHRQSLVRSFRDGTAEVETREKEPPGWGQGQGEVGNNSLGLPQGK R PASPALLPPCMLKSPKRICPVPPVSTVTGDENISDSEIEQYL RTPQEV RDFQRAQAARQAATSVPNPP</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	46.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_060780](#)

Locus ID: 55290

UniProt ID: [Q9HAW0](#)

RefSeq Size: 2007

Cytogenetics: 8p11.23


RefSeq ORF: 1257

Synonyms: BRFU; TFIIB50

Summary: This gene encodes one of the multiple subunits of the RNA polymerase III transcription factor complex required for transcription of genes with promoter elements upstream of the initiation site. The product of this gene, a TFIIB-like factor, is directly recruited to the TATA-box of polymerase III small nuclear RNA gene promoters through its interaction with the TATA-binding protein. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified BRF2 protein (Cat# TP302314). The protein was produced from HEK293T cells transfected with BRF2 cDNA clone (Cat# [RC202314]) using MegaTran 2.0 (Cat# [TT210002]).