

## Product datasheet for TP302352

### LSM14A (NM\_015578) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human LSM14A, SCD6 homolog A ( <i>S. cerevisiae</i> ) (LSM14A), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC202352 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MSGGTPYIGSKISLISKAEIRYEGILYITIDTENSTVALAKVRSFGTEDRPTDRPIPPRDEVFEYIIFRGS  
DIKDLTVCEPPKPQCCLPQDPAIVQSSLGSSTSSFQSMGSYGPFGRMPTYSQFSPSSLVGGQFGAVGVAG  
SSLTSFGTETSNSGTLPQSSAVGSAFTQDTRSLKTQLSQGRSSQLDPLRKSPTMEQAVQTASAHLPPA  
AVGRRSPVSTRPLPSASQKAGENQEHRRAEVHKVSRPENEQLRNDNKRQVAPGAPSAPRRRGRGGHRGGRG  
RFGIRRDGPMKFEKDFDFESANAQFNKEEIDREFHNKLLKEDKLEKQEKPVNGEDKGDGSDVTQNSEGN  
ADEEDPLGPNCCYDKTKSFFDNISCCDDNRERRPTWAEERRLNAETFGIPLRPNRGRGGYRGRGGLGFRGG  
RGRGGGRGGTFTAPRGRGGFRGGGRGGREFADFEYRKDNKVA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	50.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq:	<a href="#">NP_056393</a>
Locus ID:	26065
UniProt ID:	<a href="#">Q8ND56</a>
RefSeq Size:	3715
Cytogenetics:	19q13.11
RefSeq ORF:	1389
Synonyms:	C19orf13; FAM61A; RAP55; RAP55A
Summary:	Sm-like proteins were identified in a variety of organisms based on sequence homology with the Sm protein family (see SNRPD2; 601061). Sm-like proteins contain the Sm sequence motif, which consists of 2 regions separated by a linker of variable length that folds as a loop. The Sm-like proteins are thought to form a stable heteromer present in tri-snRNP particles, which are important for pre-mRNA splicing.[supplied by OMIM, Mar 2008]

### Product images:



Coomassie blue staining of purified LSM14A protein (Cat# TP302352). The protein was produced from HEK293T cells transfected with LSM14A cDNA clone (Cat# [RC202352]) using MegaTran 2.0 (Cat# [TT210002]).