

Product datasheet for **TP302378M**

Cytochrome b5 (CYB5A) (NM_148923) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cytochrome b5 type A (microsomal) (CYB5A), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202378 protein sequence Red =Cloning site Green =Tags(s)

MAEQSDEAVKYYTLEEIQKHNSKSTWLILHHKVVYDLTKFLEEHPGGEEVLRQAGGDATENFEDVGHST
DAREMSKTFIIGELHPDDRPKLNKPPETLITTIDSSSSWWTNWWIPASAVAVALMYRLYMAED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	15.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_683725
Locus ID:	1528
UniProt ID:	P00167 , A0A384ME44
RefSeq Size:	850



[View online »](#)

Cytogenetics: 18q22.3


RefSeq ORF: 402

Synonyms: CYB5; MCB5; METAG

Summary: The protein encoded by this gene is a membrane-bound cytochrome that reduces ferric hemoglobin (methemoglobin) to ferrous hemoglobin, which is required for stearyl-CoA-desaturase activity. Defects in this gene are a cause of type IV hereditary methemoglobinemia. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified CYB5A protein (Cat# [TP302378]). The protein was produced from HEK293T cells transfected with CYB5A cDNA clone (Cat# [RC202378]) using MegaTran 2.0 (Cat# [TT210002]).