

## Product datasheet for **TP302401**

### **C9orf156 (TRMO) (NM\_016481) Human Recombinant Protein**

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 9 open reading frame 156 (C9orf156), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202401 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MRGLEEPGPRPTATPCGCVKPALETGNLLTEPVGYLESCFSAKNGTPRQPSICSYSRACLRIKRIFNNP EHSMLGLEQFSHWILFVFHKNHLSCKAKVQPPRLNGAKTGVFSTRSPHRPNAIGLTLAKLEKVEGGAI YLSGIDMIHGTPVLDIKPYIAEYDSPQNVMEPLADFNLQNNQHTPNTVSQSDSKTDSCDQRQLSGCDEP Q PHHSTKRKPKCPEDRTSEENYLTHSDTARIQQAFPMHREIAVDFGLESRRDQSSSVAEEQIGPYCPEKSF SEKGTDKKLERVEGA AVLQGSRAETQPMAPHCPAGRADGAPRSVVPWVTEAPVATLEVRFTPHAEMDL G QLSSQDVGQASFKYFQSAEEAKRAIEAVLSADPRSVYRRKLCQDRLFYFTVDIAHVTCWFGDGF AEVLRI KPASEPVMHTGTPVGSLSLGS
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	48.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u><a href="#">NP_057565</a></u>
<b>Locus ID:</b>	51531
<b>UniProt ID:</b>	<u><a href="#">Q9BU70</a></u>
<b>RefSeq Size:</b>	1664
<b>Cytogenetics:</b>	9q22.33
<b>RefSeq ORF:</b>	1323
<b>Synonyms:</b>	C9orf156; HSPC219; NAP1
<b>Summary:</b>	S-adenosyl-L-methionine-dependent methyltransferase responsible for the addition of the methyl group in the formation of N6-methyl-N6-threonylcarbamoyladenosine at position 37 (m(6)t(6)A37) of the tRNA anticodon loop of tRNA(Ser)(GCU) (PubMed:25063302). The methyl group of m(6)t(6)A37 may improve the efficiency of the tRNA decoding ability. May bind to tRNA (By similarity).[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified TRMO protein (Cat# TP302401). The protein was produced from HEK293T cells transfected with TRMO cDNA clone (Cat# [RC202401]) using MegaTran 2.0 (Cat# [TT210002]).