

## Product datasheet for **TP302535**

### **C16orf48 (ENKD1) (NM\_032140) Human Recombinant Protein**


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 16 open reading frame 48 (C16orf48), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202535 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MCEGPSRISGPIPPDPTLCPDNYRRPTSAQGRLEGNAKLDLLTSDRALDTTAPRGPCIGPGAGEILERG QRGVGDVLLQLEGISLPGASLKRKDPKDHEKENLRRIREIQKRFREQERSREQGQPRPLKALWRSPKYD KVESRVKAQLQEPGPASGTESAHLRAHSRCGPGLPVSSPQPTPPGPEAKEPGLGVDFIRHNARAANK RAPRRHSCSLQVLAQVLEQQRQAQEHYNATQKGHVPHYLLERRDLWRREAARKQSQPDAMP PGHT RMP ENQRLETLTKLLQSQSLLRELVLPPAGADSLRAQSHRAELDRKLVQVEEAIKIFSRPKVFVKMDD</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	38.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_115516</a></u>
Locus ID:	84080



[View online »](#)

UniProt ID: Q9H0I2  
 RefSeq Size: 1628  
 Cytogenetics: 16q22.1  
 RefSeq ORF: 1038  
 Synonyms: C16orf48; DAKV6410; FBB11

**Product images:**


Coomassie blue staining of purified ENKD1 protein (Cat# TP302535). The protein was produced from HEK293T cells transfected with ENKD1 cDNA clone (Cat# [RC202535]) using MegaTran 2.0 (Cat# [TT210002]).