

Product datasheet for **TP302568**

RBM19 (NM_016196) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human RNA binding motif protein 19 (RBM19), 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC202568 protein sequence
Red=Cloning site Green=Tags(s)

MSRLIVKNLPLNGMKEERFRQLFAAFGLTDCSLKFTKDGKFRKFGFIGFKSEEEAQAQKHFNKSFIDTS
RITVEFCKSFGDPAKPRAWKHAQKPSQPKQPPKDDSTTPEIKKDEKKKVVAGQLEKLKEDTEFQEFSLVH
QRRQAATWANDGLDAEPSKGGKSPASDYLNFDSDSGQSEEEGAGEDLEEEASLEPKAAVQKELSDMDY
LKSKMVKAGSSSSSEEESEDEAVHCDEGSEAEEDSSATPVLQERDSRGAGQEQGMPAGKKRPPPEARAE
TEKPANQKEPTTCHTVKLRGAPFNVTEKNVMEFLAPLKPVAIRIVRNAHGNKTGYIFVDFSNEEEVKQAL
KCNREYMGGRYIEVFREKNVPTTKGAPKNTTKSWQGRILGENEEEEEDLAESGRLFVRNLPYTSTEEDLEK
LFSKYGPLSELHYPIDSLTKKPKGFATIFMFPEHAVKAYSEVDGQVFQGRMLHVLPTSTIKKEASEDASA
LGSSSYKKKKEAQDKANSASSHNWNTLFMGPNAVADAIAQKYNATKSQVFDHETKGSVAVRVALGETQLV
QEVRRFLIDNGVSLDSFSQAAAERSKTVILVKNLPAGTLAAELQETFGRFGSLGRVLLPEGGTTAIVEFL
EPLARKAFRHLAYSKFHHVPLYLEWAPVGVFSSAAPQKKKLQDTPSEPMEKDPAPETVPDGETPEDEN
PTEEGADNSSAKMEEEEEEEEEESLPGCTLFIKNLNFDTTEEKLEKVFSGVTKVSCSISKKKKNKAGV
LLSMGFVFEYRKPEQAQKALKQLQGHVVDGHKLEVRISERATKPAVTLARKKQVPRKQTTSKILVRNIP
FQAHSREIRELFSTFGELKTVRLPKKMTGTGTHRGFGVDFLTKQDAKRAFNALCHSTHLYGRRLLVLEWA
DSEVTLQALRRKTAHFHEPPKKKRSVVLDEILEQLEGSDDSDSEEQTLQL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag: C-Myc/DDK
Predicted MW: 107.2 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057280
Locus ID:	9904
UniProt ID:	Q9Y4C8
RefSeq Size:	4217
Cytogenetics:	12q24.13-q24.21
RefSeq ORF:	2880
Synonyms:	Mrd1
Summary:	This gene encodes a nucleolar protein that contains six RNA-binding motifs. The encoded protein may be involved in regulating ribosome biogenesis. Multiple alternatively spliced variants, encoding the same protein, have been identified.[provided by RefSeq, Apr 2009]

Product images:



Coomassie blue staining of purified RBM19 protein (Cat# TP302568). The protein was produced from HEK293T cells transfected with RBM19 cDNA clone (Cat# [RC202568]) using MegaTran 2.0 (Cat# [TT210002]).