

## Product datasheet for **TP302583L**

### ALKBH7 (NM\_032306) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human alkB, alkylation repair homolog 7 (E. coli) (ALKBH7), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202583 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MAGTGLLALRTLPGPSWVRGSGPSVLSRLQDAAVVRPGFLSTAEETLSRELEPELRRRRYEYDHWDAI HGFRETEKSRWSEASRAILQRVQAAAFGPGQTLLSSVHVLDEARGYIKPHVDSIKFCGATIAGLSLLSP SVMRLVHTQEPGEWLELLLEPGSLYLGRSARYDFSHEILRDEESFFGERRIPRGRISVICRSLPEGMG PGESGQPPAC
	<b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_115682</a></u>
Locus ID:	84266
UniProt ID:	<u><a href="#">Q9BT30</a></u>



[View online »](#)

RefSeq Size:	1331
Cytogenetics:	19p13.3
RefSeq ORF:	663
Synonyms:	ABH7; SPATA11; UNQ6002
Summary:	<p>May function as protein hydroxylase; can catalyze auto-hydroxylation at Leu-110 (in vitro), but this activity may be due to the absence of the true substrate (PubMed:25122757). Required to induce programmed necrosis in response to DNA damage caused by cytotoxic alkylating agents. Acts by triggering the collapse of mitochondrial membrane potential and loss of mitochondrial function that leads to energy depletion and cell death (PubMed:23666923). ALKBH7-mediated necrosis is probably required to prevent the accumulation of cells with DNA damage (PubMed:23666923). Does not display DNA demethylase activity (PubMed:23666923). Involved in fatty acid metabolism (By similarity). [UniProtKB/Swiss-Prot Function]</p>

### Product images:



Coomassie blue staining of purified ALKBH7 protein (Cat# [TP302583]). The protein was produced from HEK293T cells transfected with ALKBH7 cDNA clone (Cat# [RC202583]) using MegaTran 2.0 (Cat# [TT210002]).