

## Product datasheet for **TP302591L**

### **C16orf57 (USB1) (NM\_024598) Human Recombinant Protein**

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human chromosome 16 open reading frame 57 (C16orf57), 1 mg

**Species:** Human

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >RC202591 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MSAAPLVGYSSSGSEDESEDGMRTRPGDGSRRGQSPLPRQRFVPDSVLNMFPGTEEGPEDDSTKHGGR  
VRTFPHERGNWATHVYVPYEAKKEEFLDLLDVLLPHAQTYVPRLVRMKVFHLSLSQSWLRHHWILPFVQA  
LKARMTSFHRFFFTANQVKIYTNQEKTRTFIGLEVTSGHAQFLDLVSEVDRVMEEFNLTTFYQDPSFHLS  
LAWCVGDARLQLEGQCLQELQAIVDGFEDAEVLLRVHTEQVRCKSGNKFFSMPK

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK

**Predicted MW:** 30.1 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_078874](#)

**Locus ID:** 79650

**UniProt ID:** [Q9BQ65](#), [A0A024R6V6](#)



[View online >](#)

RefSeq Size: 2287


Cytogenetics: 16q21

RefSeq ORF: 795

Synonyms: C16orf57; hUsb1; HVSL1; Mpn1; PN

**Summary:** This gene encodes a protein with several conserved domains, however, its exact function is not known. Mutations in this gene are associated with poikiloderma with neutropenia (PN), which shows phenotypic overlap with Rothmund-Thomson syndrome (RTS) caused by mutations in the RECQL4 gene. It is believed that this gene product interacts with RECQL4 protein via SMAD4 proteins, explaining the partial clinical overlap between PN and RTS. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Mar 2011]

### Product images:



Coomassie blue staining of purified USB1 protein (Cat# [TP302591]). The protein was produced from HEK293T cells transfected with USB1 cDNA clone (Cat# [RC202591]) using MegaTran 2.0 (Cat# [TT210002]).