

## Product datasheet for **TP302606L**

### CHCHD5 (NM\_032309) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil-helix-coiled-coil-helix domain containing 5 (CHCHD5), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202606 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MQAALEVTARYCGRELEQYGQCVAAPESWQRDCHYLKMSIAQCTSSHPPIRQIRQACAQPFEAFEECLR QNEAAVGNCAEHMRRFLQCAEQVQPPRSPATVEAQPLPAS  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	12.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_115685</a></u>
Locus ID:	84269
UniProt ID:	<u><a href="#">Q9BSY4</a></u> , <u><a href="#">A0A2U9EWT9</a></u>
RefSeq Size:	698



[View online »](#)

Cytogenetics: 2q14.1  
RefSeq ORF: 330  
Synonyms: C2orf9; CHTM1; MIC14; MIX14

**Product images:**

Coomassie blue staining of purified CHCHD5 protein (Cat# [TP302606]). The protein was produced from HEK293T cells transfected with CHCHD5 cDNA clone (Cat# [RC202606]) using MegaTran 2.0 (Cat# [TT210002]).