

Product datasheet for TP302627L

PTEN (NM_000314) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--------------------------------------------------------------------------|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human phosphatase and tensin homolog (PTEN), 1 mg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC202627 protein sequence Red=Cloning site Green=Tags(s) |

MTAIIKEIVSRNKRRYQEDGFDLDTIYIPNIIAMGFPAERLEGVYRNNIDDVVRFLDSKHKNHYKIYNL
CAERHYDTAKFNCRVAQYPFEDHNPPQLELIKPFCELDQWLSDDNHVAAIHCKAGKGRGVMICAYLL
HRGKFLKAQEALDFYGEVTRDKKGVTTIPSQRRYVYYSYLLKNHLDYRPVALLFHKMMFETIPMFSGGT
CNPQFVVCQLKVKIYSSNSGPTRRDKFMYFEFPQPLPVCQDIKVEFFHKQNKMLKDKMFHFWNTFFI
PGPEETSEKVENGLCDQEIDSICSIERADNDKEYLVLTLTKNLDKANKDKANRYFSPNFKVKLYFTKT
VEEPSNPEASSSTSVTPDVSDNEPDHYRSDTTDSDPENEPFDEDQHTQITKV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


| | |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| Tag: | C-Myc/DDK |
| Predicted MW: | 47 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_000305 |
| Locus ID: | 5728 |



[View online »](#)

| | |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| UniProt ID: | P60484 , F6KD01 |
| RefSeq Size: | 5572 |
| Cytogenetics: | 10q23.31 |
| RefSeq ORF: | 1209 |
| Synonyms: | 10q23del; BZS; CWS1; DEC; GLM2; MHAM; MMAC1; PTEN1; PTENbeta; TEP1 |
| Summary: | <p>This gene was identified as a tumor suppressor that is mutated in a large number of cancers at high frequency. The protein encoded by this gene is a phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase. It contains a tensin like domain as well as a catalytic domain similar to that of the dual specificity protein tyrosine phosphatases. Unlike most of the protein tyrosine phosphatases, this protein preferentially dephosphorylates phosphoinositide substrates. It negatively regulates intracellular levels of phosphatidylinositol-3,4,5-trisphosphate in cells and functions as a tumor suppressor by negatively regulating AKT/PKB signaling pathway. The use of a non-canonical (CUG) upstream initiation site produces a longer isoform that initiates translation with a leucine, and is thought to be preferentially associated with the mitochondrial inner membrane. This longer isoform may help regulate energy metabolism in the mitochondria. A pseudogene of this gene is found on chromosome 9. Alternative splicing and the use of multiple translation start codons results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Feb 2015]</p> |
| Protein Families: | Druggable Genome, Phosphatase |
| Protein Pathways: | Endometrial cancer, Focal adhesion, Glioma, Inositol phosphate metabolism, Melanoma, p53 signaling pathway, Pathways in cancer, Phosphatidylinositol signaling system, Prostate cancer, Small cell lung cancer, Tight junction |

Product images:



Coomassie blue staining of purified PTEN protein (Cat# [TP302627]). The protein was produced from HEK293T cells transfected with PTEN cDNA clone (Cat# [RC202627]) using MegaTran 2.0 (Cat# [TT210002]).