

## Product datasheet for TP302689M

### PREI3 (MOB4) (NM\_199482) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens MOB1, Mps One Binder kinase activator-like 3 (yeast) (MOBK13), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202689 representing NM_199482 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MVMAEGTAVLRRNRPGTKAQDFYNWPDESFDMDSTLAVQQYIQQNIRADCSNIDKILEPPEGQDEGVWK  
YEHLRQFCLELNLAVKLGSECHPDTCTQMTATEQWIFLCAAHKTPKECPAIDYTRHTLDGAACLLNSNK  
YFPSRVSKESSVAKLGSVCRRIRIFSHAYFHHRQIFDEYENETFLCHRFTKFMKYNLMSKDNLIVPI  
LEEEVQNSVSGESEA

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	22.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_955776</a>
Locus ID:	25843



[View online »](#)

UniProt ID: [Q9Y3A3](#), [Q9Y3A3-2](#), [Q9UI66](#)

RefSeq Size: 2698


Cytogenetics: 2q33.1

RefSeq ORF: 675

Synonyms: 2C4D; CGI-95; MOB1; MOB3; MOBKL3; PHOCN; PREI3

**Summary:** This gene was identified based on its similarity with the mouse counterpart. Studies of the mouse counterpart suggest that the expression of this gene may be regulated during oocyte maturation and preimplantation following zygotic gene activation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. Naturally occurring read-through transcription occurs between this locus and the neighboring locus HSPE1.[provided by RefSeq, Feb 2011]

### Product images:



Coomassie blue staining of purified MOB4 protein (Cat# [TP302689]). The protein was produced from HEK293T cells transfected with MOB4 cDNA clone (Cat# [RC202689]) using MegaTran 2.0 (Cat# [TT210002]).