

Product datasheet for **TP302704L**

MSMB (NM_002443) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human microseminoprotein, beta- (MSMB), transcript variant PSP94, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202704 protein sequence Red=Cloning site Green=Tags(s)

MNVLLGSVIFATFVTLNASCYFIPNEGVPDSTRKCMDLKGKHPINSEWQTDNCETCTCYETEISCC
TLVSTPVGYDKDNCQRIFKKEDCKYIVVEKKDPKKTCSVSEWII

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	10.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002434
Locus ID:	4477
UniProt ID:	P08118
RefSeq Size:	503



[View online »](#)

Cytogenetics: 10q11.22


RefSeq ORF: 342

Synonyms: HPC13; IGBF; MSP; MSPB; PN44; PRPS; PSP; PSP-94; PSP57; PSP94

Summary: The protein encoded by this gene is a member of the immunoglobulin binding factor family. It is synthesized by the epithelial cells of the prostate gland and secreted into the seminal plasma. This protein has inhibin-like activity. It may have a role as an autocrine paracrine factor in uterine, breast and other female reproductive tissues. The expression of the encoded protein is found to be decreased in prostate cancer. Two alternatively spliced transcript variants encoding different isoforms are described for this gene. The use of alternate polyadenylation sites has been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Secreted Protein, Transmembrane

Product images:



Coomassie blue staining of purified MSMB protein (Cat# [TP302704]). The protein was produced from HEK293T cells transfected with MSMB cDNA clone (Cat# [RC202704]) using MegaTran 2.0 (Cat# [TT210002]).