

Product datasheet for **TP302715**

NDUFA2 (NM_002488) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 2, 8kDa (NDUFA2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202715 protein sequence Red =Cloning site Green =Tags(s)
	MAAAAASRGVGAKLGLREIRIHLRCRSPGSQGVRDFIEKRYVELKKANPDLPIILIRECSVQPKLWARYA FGQETNVPLNNSADQVTRALENVLSGKA
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002479
Locus ID:	4695
UniProt ID:	O43678
RefSeq Size:	726



[View online »](#)

Cytogenetics: 5q31.3


RefSeq ORF: 297

Synonyms: B8; CD14; CIB8; MC1DN13

Summary: The encoded protein is a subunit of the hydrophobic protein fraction of the NADH:ubiquinone oxidoreductase (complex 1), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane, and may be involved in regulating complex I activity or its assembly via assistance in redox processes. Mutations in this gene are associated with Leigh syndrome, an early-onset progressive neurodegenerative disorder. Alternative splicing results in multiple transcript variants.[provided by RefSeq, May 2010]

Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images:



Coomassie blue staining of purified NDUFA2 protein (Cat# TP302715). The protein was produced from HEK293T cells transfected with NDUFA2 cDNA clone (Cat# [RC202715]) using MegaTran 2.0 (Cat# [TT210002]).