

## Product datasheet for **TP302749M**

### NFYB (NM\_006166) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nuclear transcription factor Y, beta (NFYB), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202749 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MTMDGDSSTTDASQLGISADYIGGSHYVIQPHDDTEDSMNDHEDTNGSKESFREQDIYLPANVARIMKN AIPQTGKIAKDAKECVQECVSEFISFITSEASERCHQEKRKTINGEDILFAMSTLGFDSYVEPLKLYLQK FREAMKGEKGIGGAVTATDGLSEELTEEAFTNQLPAGLITTDGQQQNVMYTTSYQQISGVQQIQFS
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	22.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_006157</a>
Locus ID:	4801
UniProt ID:	<a href="#">P25208</a> , <a href="#">A0A024RBG7</a>
RefSeq Size:	3482



[View online »](#)

Cytogenetics: 12q23.3

RefSeq ORF: 621


Synonyms: CBF-A; CBF-B; HAP3; NF-YB

**Summary:** The protein encoded by this gene is one subunit of a trimeric complex, forming a highly conserved transcription factor that binds with high specificity to CCAAT motifs in the promoter regions in a variety of genes. This gene product, subunit B, forms a tight dimer with the C subunit, a prerequisite for subunit A association. The resulting trimer binds to DNA with high specificity and affinity. Subunits B and C each contain a histone-like motif. Observation of the histone nature of these subunits is supported by two types of evidence; protein sequence alignments and experiments with mutants. [provided by RefSeq, Jul 2008]

**Protein Families:** Transcription Factors

**Protein Pathways:** Antigen processing and presentation

### Product images:



Coomassie blue staining of purified NFYB protein (Cat# [TP302749]). The protein was produced from HEK293T cells transfected with NFYB cDNA clone (Cat# [RC202749]) using MegaTran 2.0 (Cat# [TT210002]).