

Product datasheet for **TP302775**

SULT1C2 (NM_001056) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sulfotransferase family, cytosolic, 1C, member 2 (SULT1C2), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202775 protein sequence Red =Cloning site Green =Tags(s)
	<p>MALTSDLGKQIKLKEVEGTLTLPATVDNWSQIQSFEAKPDDLICTYKAGTTWIQEIVDMIEQNGDVEK CQRAIIQHRHPFIEWARPPQPSGVEKAKAMPSPRILKTHLSTQLLPPSFWENNCKFLYVARNAKDCMVS YHFQRMNHMLPDPGTWEEYFETFINGKVVWGSWFDHVKGWEMKDRHQILFLFYEDIKRPKHEIRKVMQ FMGKKVDETVLDKIVQETSFEKMKENPMTNRSTVSKSILDQSISSFMRKGTVDWKNHFTVAQNERFDEI YRRKMEGTSINFCMEL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	34.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001047
Locus ID:	6819



[View online »](#)

UniProt ID:	O00338
RefSeq Size:	2799
Cytogenetics:	2q12.3
RefSeq ORF:	888
Synonyms:	humSULTC2; ST1C1; ST1C2; SULT1C1
Summary:	Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. This gene encodes a protein that belongs to the SULT1 subfamily, responsible for transferring a sulfo moiety from PAPS to phenol-containing compounds. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified SULT1C2 protein (Cat# TP302775). The protein was produced from HEK293T cells transfected with SULT1C2 cDNA clone (Cat# [RC202775]) using MegaTran 2.0 (Cat# [TT210002]).