

Product datasheet for TP302804L

Centrin 3 (CETN3) (NM_004365) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human centrin, EF-hand protein, 3 (CDC31 homolog, yeast) (CETN3), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202804 protein sequence Red=Cloning site Green=Tags(s)
	MSLALRSELVVDKTKRKKRRELSEEQKQEIKDAPFLFDTKDEAIDYHELKVAMRALGFDVKKADVLKIL KDYDREATGKITFEDFNQVTDWILERDPHEEILKAFKLFDDDDSGKISLRNLRVARELGENMSDEELR AMIEEFDKGDGGEINQEEFIAIMTGD
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004356
Locus ID:	1070
UniProt ID:	O15182



[View online »](#)

RefSeq Size: 1374

Cytogenetics: 5q14.3


RefSeq ORF: 501

Synonyms: CDC31; CEN3

Summary: The protein encoded by this gene contains four EF-hand calcium binding domains, and is a member of the centrin protein family. Centrins are evolutionarily conserved proteins similar to the CDC31 protein of *S. cerevisiae*. Yeast CDC31 is located at the centrosome of interphase and mitotic cells, where it plays a fundamental role in centrosome duplication and separation. Multiple forms of the proteins similar to the yeast centrin have been identified in human and other mammalian cells, some of which have been shown to be associated with centrosome fractions. This protein appears to be one of the most abundant centrins associated with centrosome, which suggests a similar function to its yeast counterpart. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified CETN3 protein (Cat# [TP302804]). The protein was produced from HEK293T cells transfected with CETN3 cDNA clone (Cat# [RC202804]) using MegaTran 2.0 (Cat# [TT210002]).