

Product datasheet for **TP302827M**

ING4 (NM_016162) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human inhibitor of growth family, member 4 (ING4), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202827 protein sequence Red =Cloning site Green =Tags(s)

MAAGMYLEHYLDSIENLPFELQRNFKLMRDLDQRTEDLKAIEDKLATEYMSSARSLSSSEKLALLKQIQE
AYGKCKEFGDDKVQLAMQTYEMVDKHIRRLDLDLARFEADLKEKQIESDDYDSSSSKGGKKGRTQKEKKA
ARARSKGKNSDEEAPKTAQKKLKLVRTSPEYGMPSVTFGVSHPDVLDMVPDPNEPTYCLCHQVSYGEMI
GCDNPDCSIEWFHFACVGLTTKPRGKWFCPRCSQERKKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	28.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057246
Locus ID:	51147



[View online »](#)

UniProt ID: [Q9UNL4](#)

RefSeq Size: 1458

Cytogenetics: 12p13.31


RefSeq ORF: 747

Synonyms: my036; p29ING4

Summary: This gene encodes a tumor suppressor protein that contains a PHD-finger, which is a common motif in proteins involved in chromatin remodeling. This protein can bind TP53 and EP300/p300, a component of the histone acetyl transferase complex, suggesting its involvement in the TP53-dependent regulatory pathway. Multiple alternatively spliced transcript variants have been observed that encode distinct proteins. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified ING4 protein (Cat# [TP302827]). The protein was produced from HEK293T cells transfected with ING4 cDNA clone (Cat# [RC202827]) using MegaTran 2.0 (Cat# [TT210002]).