

Product datasheet for **TP302830**

MRPS30 (NM_016640) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mitochondrial ribosomal protein S30 (MRPS30), nuclear gene encoding mitochondrial protein, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202830 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAAARCWRPLLGRPRSLHTAANAAATATETTCQDVAATPVARYPPIVASMTADSKAARLRRIERWQATV HAAESVDEKLRLTKMQFMKYMVYPQTFALNADRWYQYFTKTVFLSGLPPPPAEPEPEPEPEPALDLA ALRAVACDCLLQEHFYLRRRRRVHRYEESVISLPFLDQLVSTLVGLLSPHPALAAAALDYRCPVHFYW VRGEEIIPRGHRRGRIDDLRYQIDDKPNNQIRISKQLAEFVPLDYSVPIEPTIKCKPDKLPLFKRQYEN HIFVGSKTADPCCYGHTQFHLLPDKLRRERLLRQNCADQIEVVFRANAIASLFAWTGAQAMYQGFWSEAD VTRPFVSQAVITDGKYFSFFCYQLNLTALTTQADQNNPRKNICWGTQSKPLYETIEDNDVKGFNDVLLQ IVHLLNRPKEEKSQLEN</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	50.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.




[View online »](#)

RefSeq:	NP_057724
Locus ID:	10884
UniProt ID:	Q9NP92
RefSeq Size:	1686
Cytogenetics:	5p12
RefSeq ORF:	1317
Synonyms:	MRP-S30; PAP; PDCD9; S30mt

Summary: Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that is similar to the chicken pro-apoptotic protein p52. Transcript variants using alternative promoters or polyA sites have been mentioned in the literature but the complete description of these sequences is not available. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified MRPS30 protein (Cat# TP302830). The protein was produced from HEK293T cells transfected with MRPS30 cDNA clone (Cat# [RC202830]) using MegaTran 2.0 (Cat# [TT210002]).